

OMB MEDITECH AG

Industrie Hofacker
 CH-9606 Bütschwil · Switzerland
 Tel: +41-71-982-7201
 Fax: +41-71-982-7202
 e-Mail: sales@omb-meditech.com



► [Home](#) ► [Deutsch](#) ► [English](#) ► [Français](#) ► [L'Italiano](#)

Vapocid®-Tray System compatible to Implant/ Surgery Cassettes of various Manufacturers

- ▶ **Vapocid®-Tray System STANDARD**
- ▶ **Vapocid®-Tray System ENDO**



The Vapocid®-Tray System offers you here even more advantages:
 The maximum sterile storage period of up to 12 month in a closed cabinet offers you the optimum method for the storage of your implant- and surgery instruments and cassettes.




Renown manufacturers, such as **Institut Straumann** (Switzerland) for example recommend the Vapocid®-Tray System for this purpose!
 For the upper part of the surgery cassette use the Vapocid®-Tray 1/1 ENDO and for the lower part the Vapocid®-Tray1/1 STANDARD.

This compatibility is an additional important contribution towards increasing the hygiene standards in your practice! After treatment and cleaning of the instruments in disinfecting solution / ultrasound, prepare both trays for sterilizing in the autoclave.

The usual work process of wrapping the cassette in sterile wrapping paper is also completely eliminated!

By utilizing the Vapocid®-Tray System for this application, you can provide the necessary hygiene protection and meet with the latest requirements.



-  WORKING-TRAY
-  WASH-TRAY
-  STERILE PACKAGING

The Vapocid®-Tray System - the practical solution for the sterilization and sterile storage of your implant instruments and surgery cassettes.

► [View / Print](#)
[Adobe Acrobat File \(pdf\)](#)

◀ [back](#)