

Vapocid®-Tray System ENDO



There are many cases where it can be technically difficult or costly to meet modern hygiene requirements – but this is where the patented Vapocid®-Tray System offers a real solution,

because it can be used in three ways:


- as a washing tray (e.g. in thermo-disinfectors),
- as a normal working tray, and
- as a sterile packaging within one tray system.

Thanks to these features, the highest economic and ecological requirements can be attained

Completely variable !

- 1/1 tray carries up to 4 ENDO instrument holders
- 1/2 tray carries up to 2 ENDO instrument holders
- each instrument holder can hold 142 instruments

Instrument holders
made of plastic for
ENDO-instruments
(easy assembly)



Part No.	Description	Version	L x W X H
38935	Vapocid®-Tray System 1/1 ENDO	stainless steel	288x189x70
38937	Vapocid®-Tray System 1/1 ENDO	plastic	288x189x78
38936	Vapocid®-Tray System 1/2 ENDO	stainless steel	189x144x70
38938	Vapocid®-Tray System 1/2 ENDO	plastic	189x144x78
38941	Instrument holder ENDO 142 compl. (1 pc.)	plastic	132x83x32
38942	Instrument holder ENDO 284 compl. 2-Pack	plastic	132x83x32 /2
38944	Instrument holder ENDO 568 compl. 4-Pack	plastic	132x83x32 /4
38717	Instrument holder Standard 2-Pack	plastic flexible	-

A convincingly simple and cheap solution for ideal sterile instrument storage. Triple benefit – the system can be used in three different ways: as a working tray, a washing tray and a sterile packaging.

OMB MEDITECH AG
Industrie Hofacker
CH-9606 Bütschwil - Switzerland
Telefon +41-71-982-7201
Telefax +41-71-982-7202
URL: www.omb-meditech.com
e-Mail: sales@omb-meditech.com